

Anti-Gars Polyclonal Antibody For Wb, If/Icc, Elisa - P41250

Item Number: CM0013467



Introduction

Rabbit polyclonal antibody targeting human GARS (Glycine--tRNA ligase), validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for IVD research on protein synthesis and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	GARS Rabbit pAb
Remarks/Alias	GARS; HMN5; CMT2D; DSMAV; GlyRS; HMN5A; SMAD1; SMAJ1
Species	Human
GeneID	2617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-370 of human GARS (NP_002038.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41250
Protein Weight	83kDa
Shipping	Ice bag