

Anti-Gars Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - P41250

Item Number: CM0027504



Introduction

KO validated rabbit polyclonal antibody targeting human GARS (Glycine--tRNA ligase). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for investigating aminoacyl-tRNA synthetase biology and Charcot-Marie-Tooth disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GARS Rabbit pAb
Remarks/Alias	GARS; HMN5; CMT2D; DSMAV; GlyRS; HMN5A; SMAD1; SMAJI; RS
Species	Human
GeneID (Human)	2617
GeneID	2617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-370 of human GARS (NP_002038.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41250
Protein Weight	83kDa
Shipping	Ice bag