


Anti-Gart Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P22102

Item Number: CM0030504



Introduction

Rabbit polyclonal antibody against human GART (trifunctional purine biosynthetic protein) validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for purine biosynthesis research.

[Learn More](#)

Parameter	Specification
Product Name	GART Rabbit pAb
Remarks/Alias	AIRS; GARS; PAIS; PGFT; PRGS; GARTF; GART
Species	Human
GeneID	2618 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human GART (NP_000810.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22102
Protein Weight	108kDa
Shipping	Ice bag