


Anti-Gas7 Polyclonal Antibody For Wb, Elisa - O60861

Item Number: CM0013568



Introduction

High-quality rabbit polyclonal antibody against human GAS7, validated for WB and ELISA. Detects GAS7 in human, mouse, and rat samples, supporting neuronal differentiation research.

[Learn More](#)

Parameter	Value
Product Name	GAS7 Rabbit pAb
Remarks/Alias	GAS7; MLL/GAS7
Species	Human
GeneID	8522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-412 of human GAS7 (NP_001124303.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60861
Protein Weight	54kDa
Shipping	Ice bag