

Anti-Gasdermin B (Full Length + C Terminal) Rabbit Polyclonal Antibody For Western Blot - Q8Tax9-6

Item Number: CM0025316

Introduction

Rabbit polyclonal antibody targeting full-length and C-terminal human Gasdermin B (GSDMB, isoform 6). Validated for Western blot, suitable for pyroptosis and granzyme-mediated cell death studies.

[Learn More](#)

Parameter	Value
Product Name	Gasdermin B □ Full Length+C terminal□ Rabbit pAb
Remarks/Alias	GSDML; PP4052; GSDMB-1; PRO2521; Gasdermin B □ Full Length+C terminal□
Species	Human
Gene ID (Human)	55876
Gene ID	55876
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	Q8TAX9-6
Protein Weight	18kDa/45kDa/46kDa/47kDa
Shipping	Ice bag