


Anti-Gata1 Monoclonal Antibody For Wb, Ip, Elisa - P15976

Item Number: CM0003945



Introduction

Rabbit monoclonal antibody against human GATA1, suitable for WB, IP, and ELISA.

Immunogen: synthetic peptide within amino acids 1-100. Research use: transcriptional regulation in erythroid development.

[Learn More](#)

Parameter	Value
Product Name	GATA1 Rabbit mAb
Remarks/Alias	GF1; GF-1; NFE1; XLTT; ERYF1; NF-E1; XLANP; XLTDA; GATA-1; HAEADA; GATA1
Species	Human
Gene ID	2623
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GATA1 (NP_002040.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P15976
Protein Weight	43kDa
Shipping	Ice bag