

## Anti-Gata2 Rabbit Polyclonal Antibody For Ihc-P, Elisa - P23769

Item Number: CM0023603



### Introduction

High-quality rabbit polyclonal antibody specifically detecting human GATA2, a key endothelial transcription factor. Validated for IHC-P and ELISA, suitable for studying transcriptional regulation in endothelial cells.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | GATA2 Rabbit pAb                                                                                                          |
| Remarks/Alias    | DCML; IMD21; NFE1B; MONOMAC; GATA2                                                                                        |
| Species          | Human                                                                                                                     |
| GeneID           | 2624                                                                                                                      |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 240-340 of human GATA2(NP_116027.2). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | IHC-P, ELISA                                                                                                              |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | P23769                                                                                                                    |
| Protein Weight   | 51kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |