

Anti-Gata3 Polyclonal Antibody For Wb And Elisa - P23771

Item Number: CM0024545



Introduction

GATA3 Rabbit polyclonal antibody validated for Western blot and ELISA, reactive with human and mouse. Ideal for studying T-cell development, Th2 differentiation, and macrophage activation.

[Learn More](#)

Parameter	Value
Product Name	GATA3 Rabbit pAb
Remarks/Alias	HDR; HDRS; GATA3
Species	Human
Gene ID	2625
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GATA3 (NP_001002295.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P23771
Protein Weight	48kDa
Shipping	Ice bag