

Anti-Gata3 Rabbit Monoclonal Antibody For Wb, If, Fc - P23771

Item Number: CM0006021



Introduction

GATA3 rabbit monoclonal antibody validated for WB, IF/ICC, IF-P, FC(intra), and ELISA. Detects endogenous GATA3 in human samples, with a predicted band size of 48 kDa. Ideal for T-cell differentiation and immune response research.

[Learn More](#)

Parameter	Value
Product Name	GATA3 Rabbit mAb
Remarks/Alias	HDR; HDRS
Species	Human
GeneID (Human)	2625
GeneID	2625
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human GATA3 (NP_001002295.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, FC (intra), ELISA
Cross Reactivity	Human
SWISS	P23771
Protein Weight	48 kDa
Shipping	Ice bag