

Anti-Gata3 Rabbit Monoclonal Antibody For Wb, Ip, Elisa, Chip, Chip-Seq - P23771

Item Number: CM0002850



Introduction

Recombinant rabbit monoclonal antibody targeting GATA3, validated for WB, IP, ELISA, ChIP, and ChIP-seq in human and mouse samples. Ideal for Th2 differentiation and macrophage reprogramming research.

[Learn More](#)

Parameter	Value
Product Name	GATA3 Rabbit mAb
Remarks/Alias	HDR; HDRS; GATA3
Species	Human
GeneID (Human)	2625
GeneID	2625
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human GATA3 (P23771).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse
SWISS	P23771
Protein Weight	48kDa
Shipping	Ice bag