

Anti-Gata3 Rabbit Polymab® For Wb, If/Icc, If-P, Ihc-P, Chip, Elisa - P23772

Item Number: CM0009454



Introduction

Recombinant rabbit monoclonal antibody targeting GATA3 (UniProt P23772). Validated for WB, IF/ICC, IF-P, IHC-P, ChIP, ELISA. Detects human, mouse, and rat GATA3. Ideal for T-cell and Th2 differentiation research.

[Learn More](#)

Parameter	Value
Product Name	Gata3 Rabbit PolymAb®
Remarks/Alias	jal; Gata-3
Species	Mouse
GeneID (Human)	2625
GeneID	14462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 135-258 of mouse GATA3 (NP_032117.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23772
Protein Weight	48 kDa
Shipping	Ice bag