


Anti-Gata4 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P43694

Item Number: CM0008916



Introduction

Rabbit monoclonal antibody against human GATA4 (P43694), suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for studying cardiac development and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	GATA4 Rabbit mAb
Remarks/Alias	TOF; ASD2; VSD1; TACHD; GATA4
Species	Human
GeneID	2626
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 343-442 of human GATA4 (NP_002043.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43694
Protein Weight	45kDa
Shipping	Ice bag