

## Anti-Gata6 Polyclonal Antibody For Wb, Elisa - Q92908

Item Number: CM0024376

### Introduction

GATA6 Rabbit Polyclonal Antibody (pAb) validated for WB and ELISA. Recognizes human and mouse GATA6, a transcription factor crucial for cardiac, gastric, and skin homeostasis. Immunogen covers aa 70-160 of human GATA6. Ideal for gene regulation studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | GATA6 Rabbit pAb                                                                                                  |
| Remarks/Alias    | GATA6                                                                                                             |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 2627                                                                                                              |
| GeneID           | 2627                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 70-160 of human GATA6 (NP_005248.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q92908                                                                                                            |
| Protein Weight   | 60kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |