

Anti-Gatad2B Rabbit Polyclonal Antibody (P66-Beta) For Wb, Ihc-P, Ip, Elisa - Q8Wxi9

Item Number: CM0029572

Introduction

Rabbit polyclonal antibody against human GATAD2B (p66-beta), validated for WB, IHC-P, IP, and ELISA. Reacts with human, mouse, rat. Suitable for chromatin remodeling and transcriptional repression research.

[Learn More](#)

Parameter	Value
Product Name	GATAD2B Rabbit pAb
Remarks/Alias	p68, GANDS, MRD18, P66beta, GATAD2B
Species	Human
Gene ID (Human)	57459
Gene ID	57459
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXI9
Protein Weight	65kDa
Shipping	Ice bag