

Anti-Gatc Polyclonal Antibody For Wb, Elisa - O43716

Item Number: CM0027958

Introduction

Rabbit polyclonal antibody against human GATC, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying mitochondrial protein synthesis and tRNA transamidation.

[Learn More](#)

Parameter	Value
Product Name	GATC Rabbit pAb
Remarks/Alias	15E1.2; COXPD42; GATC
Species	Human
Gene ID (Human)	283459
Gene ID	283459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human GATC (NP_789788.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O43716
Protein Weight	15kDa
Shipping	Ice bag