


Anti-Gatm Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P50440

Item Number: CM0005608



Introduction

High-specificity anti-GATM rabbit monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat GATM (~48 kDa), a key enzyme in creatine biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	GATM Rabbit mAb
Remarks/Alias	AT; RFS; AGAT; FRTS; CCDS3; FRTS1; GATM
Species	Human
GeneID (Human)	2628
GeneID	2628
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-423 of human GATM (NP_001473.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50440
Protein Weight	48kDa
Shipping	Ice bag