


Anti-Gatm Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50440

Item Number: CM0027948



Introduction

Rabbit polyclonal antibody targeting human GATM (glycine amidinotransferase), validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for creatine biosynthesis and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	GATM Rabbit pAb
Remarks/Alias	AT, RFS, AGAT, FRTS, CCDS3, FRTS1, GATM
Species	Human
GeneID	2628
Immunogen	Recombinant protein corresponding to amino acids 38-423 of human GATM (NP_001473.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50440
Protein Weight	48kDa
Shipping	Ice bag