

Anti-Gba2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hcg7

Item Number: CM0014410

Introduction

Anti-GBA2 rabbit polyclonal antibody validated for WB and ELISA. Detects human and mouse GBA2 with a predicted molecular weight of 105 kDa. Ideal for studies on non-lysosomal glucosylceramidase function in lipid metabolism, bile acid processing, and neurodegeneration.

[Learn More](#)

Parameter	Specification
Product Name	GBA2 Rabbit pAb
Remarks/Alias	AD035, SPG46, NLGase, GBA2
Species	Human
GeneID (Human)	57704
GeneID	57704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 708-927 of human GBA2 (NP_065995.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HCG7
Protein Weight	105kDa
Shipping	Ice bag