

Anti-Gba3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H227

Item Number: CM0025999



Introduction

CamelBio supplies a rabbit polyclonal antibody against human cytosolic beta-glucosidase (GBA3). Validated for WB, IF/ICC, ELISA; cross-reacts with human and mouse. Suitable for glycosylceramide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GBA3 Rabbit pAb
Remarks/Alias	CBG; GLUC; KLRP; CBGL1; GBA3
Species	Human
GeneID	57733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-162 of human GBA3 (NP_001121904.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H227
Protein Weight	54kDa
Shipping	Ice bag