

Anti-Gbf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92538

Item Number: CM0013677

Introduction

Rabbit polyclonal antibody targeting human GBF1 (Golgi-specific brefeldin A-resistance guanine nucleotide exchange factor 1), suitable for Western blot and ELISA. Cross-reacts with mouse GBF1.

[Learn More](#)

Parameter	Value
Product Name	GBF1 Rabbit pAb
Remarks/Alias	CMT2GG; CMTDI2; CMTDIA; ARF1GEF; GBF1
Species	Human
Gene ID (Human)	8729
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human GBF1 (NP_004184.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92538
Protein Weight	206kDa
Shipping	Ice bag