


Anti-Gbp2 Polyclonal Antibody For Wb, If/Icc, Elisa - P32456

Item Number: CM0016871



Introduction

Rabbit polyclonal antibody targeting human GBP2, suitable for WB, IF/ICC, and ELISA applications. Validated in human and mouse samples, recognizing a 67 kDa protein. Ideal for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	GBP2 Rabbit pAb
Remarks/Alias	GBP2
Species	Human
Gene ID (Human)	2634
Gene ID	2634
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32456
Protein Weight	67kDa
Shipping	Ice bag