

Anti-Gbx2 Polyclonal Antibody For Wb, Elisa - P52951

Item Number: CM0014871

Introduction

Rabbit polyclonal antibody against homeobox protein GBX-2, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Suitable for studying transcription factor function in embryogenesis and development. - P52951

[Learn More](#)

Parameter	Value
Product Name	GBX2 Rabbit pAb
Remarks/Alias	GBX2
Species	Human
GeneID (Human)	2637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human GBX2 (NP_001476.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52951
Protein Weight	37kDa
Shipping	Ice bag