


## Anti-Gc1Q R/C1Qbp Polyclonal Antibody For Wb, If/Icc, Elisa - Q07021

Item Number: CM0020813



### Introduction

High-quality rabbit polyclonal antibody targeting human GC1q R/C1QBP (p32, HABP1, gC1qR). Validated for Western blot, immunofluorescence/ICC, and ELISA. Shows cross-reactivity with human, mouse, and rat. Ideal for complement system, mitochondrial function, and host-pathogen interaction studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                          |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | GC1q R/C1QBP Rabbit pAb                                                                                                                        |
| Remarks/Alias    | p32; HABP1; gC1qR; GC1QBP; SF2p32; gC1Q-R; COXPD33; SF2AP32; GC1q R/C1QBP                                                                      |
| Species          | Human                                                                                                                                          |
| GeneID (human)   | 708                                                                                                                                            |
| GeneID           | 708                                                                                                                                            |
| Immunogen        | Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 74-282 of human GC1q R/C1QBP (NP_001203.1)) |
| Source           | Rabbit                                                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                                                          |
| Application      | WB, IF/ICC, ELISA                                                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                              |
| SWISS            | Q07021                                                                                                                                         |
| Protein Weight   | 31kDa                                                                                                                                          |
| Shipping         | Ice bag                                                                                                                                        |