


Anti-Gc1Q R/C1Qbp Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q07021

Item Number: CM0003638



Introduction

Rabbit monoclonal antibody targeting human GC1q R/C1QBP (p32, HABP1). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying inflammation, complement cascade, and mitochondrial biology.

[Learn More](#)

Parameter	Value
Product Name	GC1q R/C1QBP Rabbit mAb
Remarks/Alias	p32; HABP1; gC1qR; GC1QBP; SF2p32; gC1Q-R; COXPD33; SF2AP32; GC1q R/C1QBP
Species	Human
Gene ID (Human)	708
Gene ID	708
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human GC1q R/C1QBP (Q07021)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07021
Protein Weight	31kDa
Shipping	Ice bag