

Anti-Gcat Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75600

Item Number: CM0016271

Introduction

Rabbit polyclonal antibody against human GCAT (2-amino-3-ketobutyrate coenzyme A ligase) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat. Mitochondrial target involved in glycine synthesis.

[Learn More](#)

Parameter	Value
Product Name	GCAT Rabbit pAb
Remarks/Alias	KBL; GCAT
Species	Human
Gene ID (Human)	23464
Gene ID	23464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-419 of human GCAT (NP_055106.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75600
Protein Weight	45kDa
Shipping	Ice bag