


Anti-Gcc2 / Gcc185 Rabbit Polyclonal Antibody For Wb, Ihc-P And Elisa - Q8Iwj2

Item Number: CM0018631



Introduction

Rabbit polyclonal antibody targeting human GCC2/GCC185 for Western blot, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for Golgi apparatus and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	GCC2 Rabbit pAb
Remarks/Alias	REN53; GCC185; RANBP2L4; GCC2
Species	Human
GeneID (Human)	9648
GeneID	9648
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GCC2 (NP_852118.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8Iwj2
Protein Weight	196kDa
Shipping	Ice bag