


Anti-Gcdh Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92947

Item Number: CM0028258



Introduction

A rabbit polyclonal antibody specifically recognizing GCDH, validated for WB, IHC-P, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying glutaryl-CoA dehydrogenase in amino acid catabolism and glutaric aciduria type I research.

[Learn More](#)

Parameter	Value
Product Name	GCDH Rabbit pAb
Remarks/Alias	GCD; ACAD5; GCDH
Species	Human
GeneID (Human)	2639
GeneID	2639
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 149-438 of human GCDH (NP_000150.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92947
Protein Weight	48kDa
Shipping	Ice bag