


Anti-Gckr Polyclonal Antibody For Wb, Elisa - Q14397

Item Number: CM0027929



Introduction

GCKR Rabbit pAb is a polyclonal antibody raised in rabbit, validated for Western blot (WB) and ELISA applications. Recognizes human GCKR, with cross-reactivity to mouse and rat. Useful for research on glucokinase regulation and carbohydrate metabolism.

[Learn More](#)

Parameter	Value
Product Name	GCKR Rabbit pAb
Remarks/Alias	GKRP; FGQTL5; GCKR
Species	Human
GeneID (Human)	2646
GeneID	2646
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GCKR (NP_001477.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14397
Protein Weight	69kDa
Shipping	Ice bag