

Anti-Gclc Monoclonal Antibody For Wb, If/Icc, Elisa - P48506

Item Number: CM0011787

Introduction

Anti-GCLC rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human and mouse GCLC, the catalytic subunit of glutamate-cysteine ligase, key in glutathione biosynthesis. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	GCLC Rabbit mAb
Remarks/Alias	GCL; GCS; GLCL; GLCLC; GCLC
Species	Human
GeneID (Human)	2729
GeneID	2729
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human GCLC (NP_001489.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48506
Protein Weight	73kDa
Shipping	Ice bag