

Anti-Gclc Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P48506

Item Number: CM0011167

Introduction

Reliable GCLC Rabbit monoclonal antibody for WB, IHC-P, ELISA applications; detects human and rat GCLC, the rate-limiting enzyme in glutathione biosynthesis, relevant to oxidative stress and anemia research.

[Learn More](#)

Parameter	Value
Product Name	GCLC Rabbit mAb
Remarks/Alias	GCL; GCS; GLCL; GLCLC; GCLC
Species	Human
GeneID (Human)	2729
GeneID	2729
Immunogen	Synthetic peptide corresponding to amino acids 538-637 of human GCLC (P48506)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P48506
Protein Weight	73kDa
Shipping	Ice bag