

Anti-Gclc Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48506

Item Number: CM0018705

Introduction

Rabbit polyclonal antibody targeting GCLC (glutamate-cysteine ligase catalytic subunit), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. For glutathione biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	GCLC Rabbit pAb
Remarks/Alias	GCL, GCS, GLCL, GLCLC, GCLC
Species	Human
GeneID (Human)	2729
GeneID	2729
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48506
Protein Weight	73kDa
Shipping	Ice bag