

Anti-Gclm Rabbit Monoclonal Antibody For Wb, Elisa - P48507

Item Number: CM0003890

Introduction

High-quality rabbit monoclonal antibody targeting human GCLM for WB and ELISA applications. Useful for investigating glutathione biosynthesis and oxidative stress pathways. Recognizes the 31 kDa regulatory subunit of glutamate-cysteine ligase.

[Learn More](#)

Parameter	Value
Product Name	GCLM Rabbit mAb
Remarks/Alias	GLCLR; GCLM
Species	Human
GeneID (Human)	2730
GeneID	2730
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 175-274 of human GCLM (P48507).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P48507
Protein Weight	31kDa
Shipping	Ice bag