

Anti-Gclm Rabbit Polyclonal Antibody For Wb And Elisa - P48507

Item Number: CM0016602

Introduction

High-quality rabbit polyclonal antibody targeting human GCLM (glutamate-cysteine ligase modifier subunit). Validated for WB and ELISA, detects GCLM in human and mouse samples. Essential tool for glutathione biosynthesis and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	GCLM Rabbit pAb
Remarks/Alias	GLCLR; GCLM
Species	Human
GeneID	2730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human GCLM (NP_002052.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48507
Protein Weight	31kDa
Shipping	Ice bag