

## Anti-Gcm1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Np62

Item Number: CM0024814

### Introduction

Rabbit polyclonal antibody targeting human GCM1, a transcription factor essential for placental development. Suitable for WB, IF/ICC, ELISA. Cross-reacts with human. Ideal for trophoblast differentiation studies.

[Learn More](#)

| Parameter        | Specification                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | GCM1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | GCMa; hGCMa; GCM1                                                                                                  |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 8521                                                                                                               |
| Gene ID          | 8521                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 207-436 of human GCM1 (NP_003634.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human                                                                                                              |
| SWISS            | Q9NP62                                                                                                             |
| Protein Weight   | 49kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |