

Anti-Gcn1L1 Polyclonal Antibody (Ko Validated) For Wb And Elisa - Q92616

Item Number: CM0023554

Introduction

Rabbit polyclonal antibody targeting human GCN1L1 (GCN1L1, PRIC295), validated for Western blot and ELISA. Ideal for ribosome collision and integrated stress response research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GCN1L1 Rabbit pAb
Remarks/Alias	GCN1L; GCN1L1; PRIC295; L1
Species	Human
GeneID (Human)	10985
GeneID	10985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2462-2671 of human GCN1L1 (NP_006827.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92616
Protein Weight	293kDa
Shipping	Ice bag