

Anti-Gcnt3 Rabbit Polyclonal Antibody For Wb And Elisa - O95395

Item Number: CM0015864

Introduction

Rabbit polyclonal antibody against human GCNT3, a mucin-type core 2/4 O-glycan branching enzyme. Validated for Western blot and ELISA. Ideal for glycobiology and mucin biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	GCNT3 Rabbit pAb
Remarks/Alias	GNTM; C24GNT; C2GNT2; C2GNTM; C2/4GnT; GCNT3
Species	Human
Gene ID (Human)	9245
Gene ID	9245
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-180 of human GCNT3 (NP_004742.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95395
Protein Weight	51kDa
Shipping	Ice bag