


Anti-Gcsh Polyclonal Antibody For Wb, If/Icc, Elisa - P23434

Item Number: CM0018157



Introduction

Rabbit polyclonal antibody against human GCSH, validated for WB, IF/ICC, and ELISA. Detects endogenous GCSH in human, mouse, and rat samples. Suitable for glycine metabolism and lipoylation studies.

[Learn More](#)

Parameter	Value
Product Name	GCSH Rabbit pAb
Remarks/Alias	GCE; NKH; GCSH
Species	Human
GeneID (Human)	2653
GeneID	2653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human GCSH (NP_004474.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23434
Protein Weight	19kDa
Shipping	Ice bag