


Anti-Gcsh Rabbit Monoclonal Antibody (Wb, If/Icc, Elisa) - P23434

Item Number: CM0009170



Introduction

High-specificity GCSH Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Detects human GCSH (P23434), a mitochondrial glycine cleavage system H protein involved in glycine metabolism and lipoylation. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	GCSH Rabbit mAb
Remarks/Alias	GCE; NKH; GCSH
Species	Human
GeneID	2653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-173 of human GCSH (NP_004474.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P23434
Protein Weight	19kDa
Shipping	Ice bag