


Anti-Gdap1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tb36

Item Number: CM0027920



Introduction

Rabbit polyclonal antibody targeting human GDAP1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for mitochondrial fission and Charcot-Marie-Tooth disease studies.

[Learn More](#)

Parameter	Value
Product Name	GDAP1 Rabbit pAb
Remarks/Alias	CMT4; CMT4A; CMTRIA; GDAP1
Species	Human
Gene ID (Human)	54332
Gene ID	54332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human GDAP1 (NP_001035808.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TB36
Protein Weight	41kDa
Shipping	Ice bag