


Anti-Gdf-8/Myostatin Rabbit Monoclonal Antibody For Wb And Elisa - O14793

Item Number: CM0005733



Introduction

Rabbit monoclonal antibody targeting human GDF-8/Myostatin (O14793), validated for Western blot (WB) and ELISA with cross-reactivity in human, mouse, and rat. Ideal for skeletal muscle growth research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	GDF-8/Myostatin Rabbit mAb
Remarks/Alias	GDF8; MSLHP; GDF-8/Myostatin
Species	Human
Gene ID (Human)	2660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-375 of human GDF-8/Myostatin (NP_005250.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14793
Protein Weight	43kDa
Shipping	Ice bag