

Anti-Gdf15 Polyclonal Antibody For If/Icc, Elisa - Q99988

Item Number: CM0021676



Introduction

This anti-GDF15 rabbit polyclonal antibody is validated for IF/ICC and ELISA, specific to human GDF15, a hormone involved in appetite regulation, cancer cachexia, and pregnancy nausea. Suitable for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	GDF15 Rabbit pAb
Remarks/Alias	PDF; MIC1; PLAB; MIC-1; NAG-1; PTGFB; GDF-15; GDF15
Species	Human
GeneID (Human)	9518
GeneID	9518
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 197-308 of human GDF15 (NP_004855.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q99988
Protein Weight	34kDa
Shipping	Ice bag