

Anti-Gdf15 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q99988

Item Number: CM0018614



Introduction

GDF15 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat GDF15/MIC-1, a stress-response hormone and cachexia marker. Ideal for cancer, metabolic, and pregnancy research.

[Learn More](#)

Parameter	Value
Product Name	GDF15 Rabbit pAb
Remarks/Alias	PDF; MIC1; PLAB; MIC-1; NAG-1; PTGFB; GDF-15; GDF15
Species	Human
GeneID	9518
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-308 of human GDF15 (NP_004855.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99988
Protein Weight	34kDa
Shipping	Ice bag