

Anti-Gdf9 Polyclonal Antibody For Wb, Elisa - O60383

Item Number: CM0026318

Introduction

Rabbit polyclonal antibody targeting human growth/differentiation factor 9 (GDF9) for Western blot and ELISA. Reacts with human, mouse, and rat. Useful in ovarian folliculogenesis research. UniProt O60383.

[Learn More](#)

Parameter	Value
Product Name	GDF9 Rabbit pAb
Remarks/Alias	POF14; GDF9
Species	Human
GeneID (Human)	2661
GeneID	2661
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 205-454 of human GDF9 (NP_005251.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60383
Protein Weight	51kDa
Shipping	Ice bag