

Anti-Gdi1 Polyclonal Antibody For Wb And Elisa - P31150

Item Number: CM0025152

Introduction

Rabbit polyclonal antibody against human GDI1 (Rab GDP dissociation inhibitor alpha), validated for WB and ELISA. Cross-reacts with human and mouse GDI1. Ideal for studying Rab protein regulation, membrane trafficking, and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	GDI1 Rabbit pAb
Remarks/Alias	1A; GDIL; MRX41; MRX48; OPHN2; XAP-4; XLID41; RABGD1A; RABGDIA; GDI1
Species	Human
GeneID (Human)	2664
GeneID	2664
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 188-447 of human GDI1 (NP_001484.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31150
Protein Weight	51kDa
Shipping	Ice bag