

Anti-Gdi2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P50395

Item Number: CM0028654

Introduction

Rabbit polyclonal antibody against human GDI2 (Rab GDP dissociation inhibitor beta) validated in WB, IF/ICC, ELISA. Suitable for human, mouse, rat samples. Offers reliable detection for intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GDI2 Rabbit pAb
Remarks/Alias	RABGDIB; HEL-S-46e; I2
Species	Human
Gene ID (Human)	2665
Gene ID	2665
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 296-445 of human GDI2 (NP_001485.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50395
Protein Weight	51kDa
Shipping	Ice bag