

Anti-Gelsolin Monoclonal Antibody For Wb, Ihc-P, Elisa - P06396

Item Number: CM0006017



Introduction

Rabbit monoclonal antibody against human Gelsolin (P06396). Validated for WB, IHC-P, ELISA. Reacts with human and mouse. Ideal for actin dynamics, cytoskeleton, and ciliogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	Gelsolin Rabbit mAb
Remarks/Alias	ADF; AGEL; Gelsolin
Species	Human
GeneID (human)	2934
GeneID	2934
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human Gelsolin (NP_000168.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06396
Protein Weight	86kDa
Shipping	Ice bag