

Anti-Gelsolin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06396

Item Number: CM0019479



Introduction

Rabbit polyclonal antibody targeting human Gelsolin (ADF, AGEL), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for actin cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	Gelsolin Rabbit pAb
Remarks/Alias	ADF; AGEL; Gelsolin
Species	Human
GeneID	2934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 282-438 of human Gelsolin (NP_000168.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06396
Protein Weight	86kDa
Shipping	Ice bag