

Anti-Gelsolin Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ihc-P, Elisa - P06396

Item Number: CM0004861

Introduction

[KD Validated] Recombinant anti-Gelsolin rabbit monoclonal antibody for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Recognizes human Gelsolin (P06396), 86 kDa. Suitable for actin cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Gelsolin Rabbit mAb
Remarks/Alias	ADF; AGEL; [KD Validated] Gelsolin
Species	Human
GeneID	2934
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-782 of human Gelsolin (NP_000168.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P06396
Protein Weight	86kDa
Shipping	Ice bag