

Anti-Gelsolin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06396

Item Number: CM0023407



Introduction

KD validated rabbit polyclonal antibody targeting human Gelsolin (ADF, AGEL). Suitable for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Gelsolin modulates actin filament dynamics.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Gelsolin Rabbit pAb
Remarks/Alias	ADF; AGEL; [KD Validated] Gelsolin
Species	Human
GeneID (Human)	2934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-782 of human Gelsolin (NP_000168.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06396
Protein Weight	86kDa
Shipping	Ice bag