

Anti-Gem Rabbit Polyclonal Antibody For Wb, Elisa - P55040

Item Number: CM0026538



Introduction

High-quality rabbit polyclonal antibody against human GTP-binding protein GEM (KIR) for Western blot and ELISA. Reacts with human samples; ideal for studying GEM-mediated signal transduction at the plasma membrane.

[Learn More](#)

Parameter	Specification
Product Name	GEM Rabbit pAb
Remarks/Alias	KIR; GEM
Species	Human
Gene ID (Human)	2669
Gene ID	2669
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GEM (NP_005252.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55040
Protein Weight	34kDa
Shipping	Ice bag